

Cytometer Configuration								
Software Version	BD FACSDiva 9.0							
Export Time	11/09/20 12:00:18 PM							
User								
Cytometer	LSRFortessa							
Serial Number	R658226R1128							
Configuration Name	X20- 2-Blue 3-Red 6-Violet 4-Yel Gr 3-UV (658226R1)							
Laser Name	Type	Laser	Power	Detector Array	Detector	Mirror	Filter	Parameter
BLUE	Custom	488	60	Octagon	A	655 LP	670/30 BP	7-AAD/PerCP/PerCP-Cy5.5
						524 LP	542/27BP	Or YFP
					B	505 LP	530/30 BP	FITC
						495 LP	510/20 BP	Or GFP
VIOLET	Custom	405	50	Octagon	C		488/10 BP	SSC
					A	750 LP	780/60 BP	BV785
					B	685 LP	710/50 BP	BV711
					C	635 LP	670/30 BP	BV650
					D	600 LP	610/20 BP	BV605
					E	505 LP	525/50 BP	BV480 or BV510
RED	Custom	640	40	Octagon	A	750 LP	780/60 BP	APC-Cy7/APC-Alexa 750
					B	690 LP	730/45 BP	APC-Alexa 700/Alex 700
					C		670/30 BP	APC
YELLGRN	Custom	561	50	Octagon	A	750 LP	780/60 BP	PE-Cy7
					B	635 LP	670/30 BP	PE-Cy5
					C	600 LP	610/20 BP	PE-Texas Red/PI
					D		585/15 BP	PE
UV	Custom	355	15	Trigon	A	685 LP	740/35 BP	BD Horizon BUV737
					B	450LP	515/30 BP	BUV496/Indo 1
						or 420LP	450/50 BP	Dapi/Hoest/FVS440
					C		379/28 BP	Brilliant UV 395